

RSAConference2016

San Francisco | February 29 – March 4 | Moscone Center

SESSION ID: MBS-W05

Wireless Infusion Pumps: Securing Hospitals' Most Ubiquitous Medical Device



#RSAC



Connect to
Protect

Nate Lesser

Deputy Directory
National Cybersecurity Center of
Excellence
NIST



POPULAR SCIENCE

HEALTH

HACKED MEDICAL DEVICES MAY BE THE BIGGEST CYBER SECURITY THREAT IN 2016

Healthcare IT News

TOPICS SIGN UP MAIN MENU

Threat matrix: Malware and hacking pose dangers to medical devices

CSO

NEWS

Attackers targeting medical devices to bypass hospital security



HealthData
Management

Lax medical device security needs urgent action

SECURITYWEEK
INTERNET AND ENTERPRISE SECURITY NEWS, INSIGHTS & ANALYSIS

Medical Devices Used as Pivot Point in Hospital Attacks: Report

Risks, Threats, Vulnerabilities to Infusion Pumps



#RSAC

RISKS

- Patient safety (lives)
- Operational/downtime
- Patient trust & staff morale

THREATS

- Collateral damage
- Malware
- Theft/loss
- Lateral attack
- Hacktivism

VULNERABILITIES

- Tightly regulated "turn-key" systems
- Long useful life
- Poorly protected & patched
- No detection & alerting
- System complexity

Infusion Pumps



- Infusion pumps are medical devices that deliver fluids and medications into a patient's body in controlled amounts.
- Infusion pumps among the most are ubiquitous connected medical devices in hospitals.
 - Estimated that more than 2 million infusion pumps are in use in hospitals.



Infusion Pumps in the News



#RSAC

WIRED

HACKER CAN SEND FATAL DOSE TO HOSPITAL DRUG PUMPS



BloombergBusiness

Hospital Drug Pump Can Be Hacked Through Network, FDA Warns

WIRED

VIDEO SHOWS A TERRIFYING DRUG INFUSION PUMP HACK IN ACTION

Infusion Pump Vulnerability Notices



#RSAC



- FDA Safety Communication on Vulnerabilities of Infusion Pump Systems, May 13, 2015: recommendation for health care facilities to secure affected infusion pumps.



INDUSTRIAL CONTROL SYSTEMS
CYBER EMERGENCY RESPONSE TEAM


- ICS-CERT Advisory (ICSA1512501B): provides details of vulnerabilities including CVE numbers and CVSS scores.
- Postmarket Management of Cybersecurity in Medical Devices: draft guidance on cybersecurity risk management, remediation, and reporting



Infusion Pump Physical Architecture



#RSAC





Demo




Exploit #1: Compromise Patient Information



#RSAC

- TELNET access to “root” account
 - No authentication required
 - CVE-2015-3459
- Risk
 - Gain root access to pump
 - Pump used as a pivot




Exploit #1: Compromise Patient Information



#RSAC

- Unauthenticated telnet to pump
 - Use pump to pivot
- Extract desired information
- Move information to web server
- Extract via web browser




Exploit #2: Crash Communication Subsystem



#RSAC

- Denial of service attack
 - Stops network communication
 - Using resource consumption
 - Can corrupt flash file system
- Risk
 - Stops pump from sending alerts and messages




Exploit #2: Crash Communication Subsystem



#RSAC

- Normal Operations
- Error condition
 - Empty vial
 - Alarm sent
- Error condition cleared
 - Nurse replaces vial
- Infusion resumes




Exploit #2: Crash Communication Subsystem



#RSAC

- Anonymous FTP to pump
- Upload gzip'd file
- TELNET access to “root” account
- gunzip file
 - Consumes all available RAM
 - Corrupts part of file system
 - No file system repair tools



All Devices Have Vulnerabilities



#RSAC

- Although the vulnerabilities might not be known
- Just because it's vulnerable doesn't mean it's exploitable
- Patching may not be possible right away



Wireless Infusion Pump Problem Statement



#RSAC

RESEARCH

- Help health care delivery organizations understand risks & secure medical devices on an enterprise network
- Focus on wireless infusion pumps as archetype of medical device

BUILD

- Assess risk
- Identify mitigating security technologies
- Build example implementation
- Independent validation of implementation

SHARE

- Publish practice guide: NIST Special Publication 1800 series

Challenges



#RSAC

- Firmware version control
 - Multiple versions in service
- Access control
 - Physical and networked
 - “Break the glass”
 - Malware or other unexpected software on pump
- Wireless access point and network configuration
- Alarms
- Asset management and monitoring
- Identity management and Credentialing
- Maintenance and firmware updates
- Pump variability
 - Multiple types of pumps
 - Multiple models in usage





MISSION

ACCELERATE ADOPTION OF
SECURE TECHNOLOGIES

Collaborate with innovators to provide
real-world, standards-based
cybersecurity capabilities that address
business needs



VISION

ADVANCE CYBERSECURITY

A secure cyber infrastructure that
inspires technological innovation
and fosters economic growth



GOAL 1

PROVIDE PRACTICAL
CYBERSECURITY

Help people secure their data and
digital infrastructure by equipping
them with practical ways to implement
standards-based cybersecurity
solutions that are modular, repeatable
and scalable



GOAL 2

INCREASE RATE OF
ADOPTION

Enable companies to rapidly deploy
commercially available cybersecurity
technologies by reducing
technological, educational and
economic barriers to adoption



GOAL 3

ACCELERATE
INNOVATION

Empower innovators to
creatively address
businesses' most pressing
cybersecurity challenges in a
state-of-the-art, collaborative
environment



Standards-based: Apply relevant local, national and international standards to each security implementation and account for each sector's individual needs; demonstrate reference designs for new standards



Modular: Develop reference designs with individual components that can be easily substituted with alternates that offer equivalent input-output specifications



Repeatable: Enable anyone to recreate the NCCoE builds and achieve the same results by providing a complete practice guide including a reference design, bill of materials, configuration files, relevant code, diagrams, tutorials and instructions



Commercially available: Work with the technology community to identify commercially available products that can be brought together in reference designs to address challenges identified by industry



Usable: Design usable blueprints that end users can easily and cost-effectively adopt and integrate into their businesses without disrupting day-to-day operations




Open and transparent: Use open and transparent processes to complete work, and seek and incorporate public comments on NCCoE documentation, artifacts and results

Participate



#RSAC




301-975-0200



hit_nccoe@nist.gov



<http://nccoe.nist.gov>



100 Bureau Dr, M/S 2002
Gaithersburg, MD 20899